


Product: Patella Tracer PF Knee Brace
Models: 64X

- Features textured C-shaped buttress for maximum patella control, versatile strapping system that provides dynamic treatment and hinges for medial/lateral support
- Indicated for patella subluxation, patello femoral syndrome and as part of a rehabilitation for patella tendon injuries where medial/lateral support is required
- Manufactured from premium grade latex-free neoprene, a neoprene alternative; compression molded textured EVA buttress and high density plastic hinges



Order #	Size	Knee Circumference
641	X-Small	12" - 13"
643	Small	13" - 14"
645	Medium	14" - 15"
647	Large	15" - 16"
648	X-Large	16" - 18"
649	2X-Large	18" - 20"
649-9	3X-Large	20" - 22"
649-99	4X-Large	22" - 24"

SIZING: Measure around the center of the knee with leg extended.

HOW TO APPLY:

Loosen straps, making sure they remain threaded through d-rings. After determining side of kneecap (patella) that requires support, position buttress on affected side and slide sleeve over knee. Make sure seam is placed posteriorly in back of leg. Your kneecap (patella) should be centered in the circular opening with buttress placed on rim of the kneecap (patella). Grasp both straps by hook tip and pull back onto loop portion of strap. Adjust straps to desired tension or comfort level.

CARE INSTRUCTIONS:

Hand wash in cold water with mild soap. Wrap in a towel to remove excess water (do not twist or wring). Air dry away from heat. Do not tumble dry. Do not bleach.

CAUTION:

Wear as recommended. Always consult a physician when encountering persistent pain, discomfort, inflammation, or swelling.

This QR code can be scanned by your smart phone and directs you to the

Product Detail Page.

SCAN

