



## Min-Knee Hinged Knee Brace™

Suggested  
HCPC L1832  
and L1833

### Features & Benefits

- Wrap around design allows for exact fit & comfort
- Provides for mild to moderate medial and lateral support in addition to the support offered by the strapping.
- Polycentric hinge provides flexion / extension.
- Stops at 0°, 15°, 30°, 60° and 90°
- Shorter knee brace design is less bulky and therefore more patient-friendly.

### Indications

Indicated for mild to moderate ligament joint pain, sprain, strains of the knee ligaments. Provides support, stabilization and compression of the knee joint following ligament or meniscus injuries.

### Sizes

Order #	Size	Knee Circumference
623	Small	12½" - 14½"
625	Medium	14" - 16"
627	Large	15½" - 18"
628	X Large	17½" - 20"
629	2X Large	19½" - 22"
629-9	3X Large	21½" - 24"

Measure around the center of the knee with leg extended

