



## The Weave™ Series

### Features & Benefits

- Single pull design
- Rigid panels provide superior posterior and Anterior support
- Patient friendly especially for arthritic patients
- Quick pull 6:1 mechanical advantage compression system
- Breathable mesh belt
- Uniform Compression
- Lightweight
- Ability to wear under cloths
- Aerospace fabric-moisture wicking and breathable material

### Indications

Spinal Stenosis, Spondylolisthesis, Spondylolysis, Chronic Back Pain, Spinal Instability, Post-Op Disk Surgery, Facet Syndrome, Degenerative, Herniated, Bulging Disc, Post-Operative Lumbar Fusion



**Weave™ 27**  
Suggested HCPC  
L0627 and L0642  
10" high back panel



**Weave™ 31**  
Suggested HCPC  
L0631 and L0648  
14" high back panel



**Weave™ 77**  
Suggested HCPC  
L0637 and L0650  
14" high back panel plus side panels for extra lateral stability

### Sizes

Weave™ 27 Order #	Weave™ 31 Order #	Weave™ 77 Order #	Size	Waist Circumference
733	743	773	Small	28" - 32"
735	745	775	Medium	32" - 36"
737	747	777	Large	36" - 40"
738	748	778	X Large	40" - 46"
739	749	779	2X Large	46" - 50"
739-9	749-9	779-9	3X Large	50" - 56"
739-99	749-99	779-99	4X Large	56" - 62"