



AO Enforcer ACL Knee Brace

Suggested
HCPC
L1845/L1852

Features & Benefits

- Designed to offer excellent support and protection for ACL, MCL and other knee ligament repairs and/or instabilities
- Light weight, high strength Aluminum Alloy frames
- State-of-the-Art low profile ICR hinge closely mimics natural knee motion
- Customizable non-elastic straps and comfort fit condylar pads for medial-lateral stability and a superior & improved fit
- Adjustable straps and a range of pads allow for fit customization
- Molded silicone pads offer non-slip suspension and avoid tibial pressure
- Variety of extension stops to limit range of motion
- Aluminum frames allows for additional customization by medical professionals



Indications

Provides stability and support following injury or reconstruction of the ACL, PCL, MCL, LCL and menisci. Can be used prophylactically to prevent hyperextension injuries.

Sizes

Order # Right	Order # Left	Size	Thigh Circumference
941-R	951-L	X-Small	13½" - 15½"
943-R	953-L	Small	15½" - 18½"
945-R	955-L	Medium	18½" - 21"
947-R	957-L	Large	21" - 23½"
948-R	958-L	X Large	23½" - 26½"
949-R	959-L	2X Large	26½" - 29"

Note: Thigh circumference should be measured approximately 6" above the knee.

Advanced Orthopaedics, Inc.
326 Main Street
Harleysville, PA 19438
800.270.7074
info@advanced-orthopaedics.com

This QR code can be scanned by your smart phone and directs you to the Product Detail Page.

