



Artic Hip Brace - Available Now!

Suggested
HCPC L1686

Features & Benefits

- Lightweight and low-profile hinge
- Works with recycling cold units while controlling hip instability
- Easily adjustable (7 points) with allen and thumb screws
- Universal/Ambidextrous – Works for both sides of the hip
- Controls ROM at 15° Increments
- Adjustable in length to patients 5'2" - 6'6"
- Customizable to various leg shares and sizes
- Convenient cold therapy pocket and straps work with or without brace
- Accommodates both anterior & posterior hip procedures



Indications

Hip replacement (arthroplasty), Hip revision surgery, Hip dislocation, Hip dysplasia, Post-op proximal hamstring repair, Repairs of labral structure with or without gluteus involvement, Pre or Post-op

Sizes

Order #	Size
ABH001	Universal Fit

