



## COOLFORCE KNEE BRACE

### Features & Benefits

- Liner and hand pump to force air into the brace and increase compression around the knee
- Removable ice pack designed to speed up recovery times and reduce pain
- Medial and lateral hinge stays provide excellent support
- 4-way stretch nylon fabric Piping along thigh and calf areas of support to stop flexing Tapered in back for comfortable fit
- Hyperextension stops at 170°
- Removable hinges for custom fitting



### Indications

Indicated for mild to moderate ligament joint pain, sprain, strains of the knee ligaments. Provides support, stabilization and compression of the knee joint following ligament or meniscus injuries.

#### Order #

5530 (w/o Hinges)  
 5535 (w/ROM Hinges)  
 5545 (w/Poly Hinge)  
 33

#### Size

Universal  
 Universal  
 Universal  
 Add'l Cold Pack

#### Sugg. HCPC

N/A  
 L1847/L1848  
 L1810/L1812  
 N/A



## THE COOLFORCE THERAPY LINE

This QR code can be scanned by your smart phone and directs you to the Product Detail Page.

