advanced-orthopaedics.com

Suggested HCPC

L1832 and L1833



Cross-Fit Universal Hinged Knee Brace™

Features & Benefits

- . Medial and lateral hinge stays provide excellent support
- Enclosed tubular donut buttress
- . Neoprene coated with 4-way stretch nylon fabric
- . Piping along thigh and calf areas of support to stop flexing
- . Tapered in back for comfortable fit
- . Hyperextension stops at 170°
- . Easy removable hinges for custom fitting
- . Retail packaging available

Indications

Indicated for mild to moderate ligament joint pain, sprain, strains of the knee ligaments. Provides support, stabilization and compression of the knee joint following ligament or meniscus injuries.

Sizes

Order # 6100

Size Universal



Retail Packaging Available

Advanced Orthopaedics, Inc. 326 Main Street Harleysville, PA 19438 800.270.7074 info@advanced-orthopaedics.com



